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NG-Test® CARBA 5

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- Yohei D et al. Carbapenemase-Producing Enterobacteriaceae. Semin Respir Crit Care Med. 2015; 36(1): 74-84.
- Nordmann P et al. Rapid Detection of Carbapenemase-producing Enterobacteriaceae. Emerging Infectious Diseases. 2012; 18(9).
- Bush K et al. Updated Functional Classification of β-Lactamases. Antimicrobial Agents and Chemotherapy. 2010; 54(3):969-976.
- Nordmann P et al. Carbapenemases in enterobacteriaceae. Archives de Pédiatrie. 2010;17:S154-S162.
- Demir Y et al. Investigation of VIM, IMP, NDM-1, KPC AND OXA-48 enzymes in Enterobacteriaceae strains. Pak J Pharm Sci. 2015; May;28(3 Suppl):1127-33.
- Oteo J et al. Evolution of carbapenemase-producing Enterobacteriaceae at the global and national level: what should be expected in the future? Enferm Infecc Microbiol Clin. 2014; 32 Suppl 4:17-23.
- Kulkarni MV et al. Use of imipenem to detect KPC, NDM, OXA, IMP, and VIM carbapenemase activity from gram-negative rods in 75 minutes using liquid chromatography-tandem mass spectrometry. J Clin Microbiol. 2014; 52(7):2500-5.
- Boutal H et al. 2018. A multiplex lateral flow immunoassay for the rapid identification of NDM-, KPC-, IMP- and VIM-type and OXA-48-like carbapenemase-producing Enterobacteriaceae. J Antimicrob Chemother.
- * Baeza LL, et al. A comparison of five methods for detection of carbapenemases in Enterobacteriaceae with proposal of a new algorithm. Clinical Microbiology and Infection.
- Hopkins KL, Meunier D, Naas T, Volland H, Woodford N. 2018. Evaluation of the NG-Test® CARBA-5 multiplex immunochromatographic assay for the detection of KPC, OXA-48-like, NDM, VIM and IMP carbapenemases. J Antimicrob Chemother 73:3523-3526
- Volland H et al., Improvement of the immunochromatographic NG-Test® CARBA-5 assay for the detection of IMP-variants previously undetected. Antimicrobial Agents and Chemotherapy. 2019
- Klavins A et al. Evaluation of NG-Test® CARBA-5 with Carbapenemase producing organisms and non-carbapenemase producing organisms. CPHM 991. ASM Microbe 2019- Chicago
- * Hazim O. et al. Comparative Evaluation of Five Assays for Detection of Carbapenemases with a Proposed Scheme for Their Precise Application, The Journal of Molecular Diagnostics 2020;22(9)

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